SINGLE
CELL
MODEL

PermeH₂O

Water Vapour

FOR FILMS

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Water vapour permeability tester with embedded controls of temperature and relative humidity.

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PermeH₂O is an instrument for testing the water vapour permeability (WVTR) through plastic films, monolayer or multilayer barrier films, metallised or with surface coating, laminated or co-extruded used especially for food, beverage, pharmaceutical and electronics packaging applications.

The dynamic-differential technique of the analysis used in PermeH₂O, avoids completely any possible erroneous measurements due to lateral leakages caused by possible surface defects or roughness of the sample. This is not possible with other techniques of analysis.

The software PermeH₂O ExtraSolution® allows the user to work with full autonomy after only few hours of training. In fact, the principal parameters such as temperature, relative humidity, and conditioning time are requested at the start of the software and then it is all automatic. In this way any human action is reduced at the minimum and any possible error caused by shortage of experience is avoided.

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
PermeH₂O complies with the norms ASTM F1249, TAPPI T557, JIS K-1729 and ISO 15106-2 for WVTR determination through plastic film and sheeting using an infrared sensor.

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The software stabilizes the parameters with high accuracy and precision over all the measurement.

Testing conditions (temperature, relative humidity, concentration,...) are shown both as numeric data and graphical representation and saved in HTML file and in text file to allow the full compatibility with any other software. All the gas flows are electronically controlled so that the instrument is insensible to changes of pressure both in the gas line and atmospheric. The single point closure system and the special design of the chamber eliminate the need of cutting the samples to be tested.

Extrasolution made by PermTech

Via del Bagno 16 55036 Pieve Fosciana (LU), Italy info@permtech.it www.permtech.it
Tel. +3905831646078 Tel. +3903311976173 Mob.+393389028861 Fax +3905831641255